v 2.6 2019

Models:

OPTIKA PRO VIEW, OPTIKA LITE VIEW

Type:

SOFTWARE

Instructions:

Before proceeding with the SW installation, please check the table below "Software Function list" to identify the most suitable software.

SOFTWARE FUNCTION LIST

	FU	NCTION	OPTIKA PRO VIEW	OPTIKA LITE VIEW
	Simultaneous management of several cameras		x	x
	GUI (Graphical User Interface)		x	
	Report generator		X	
	Archiving		X	х
	Language	Catalan	X	x
		Chinese (simpl.)	X	х
		Chinese (trad.)	X	х
		Korean	X	х
A A		English	X	X
R		French	X	x
GENERAL		German	X	X
8		Indonesian	X	х
		Italian	X	х
		Japanese	X	X
		Polish	х	x
		Russian	х	х
		Spanish	X	X
		Swedish		
		Thai	X	Х
		Turkish	x	Х

	FU	NCTION	OPTIKA PRO VIEW	OPTIKA LITE VIEW
ENTS	Measurements on "live"		x	
	Measurements on "captured"		X	
	2D Measurements	Line	х	
		Angle	x	
		Parallel lines	X	
Ш		Rectangle	x	
MEASUREMI		Ellipse	х	
		Circle	x	
		Annulus	X	
		Arc	x	
		Curve	X	
		Polygon	x	
	Particle count		х	
	Export to Excel		X	

4083.Wifi, **4083.4** and **4083.EC2** work with Vision Lite only.

C-E2 work with PROVIEW and LITEVIEW only.

Cameras with HDMI connection only, do not require any software; C-HP works with PROVIEW only, isn't compatible with Linux and MAC OS.



SOFTWARE FUNCTION LIST

FUN	ICTION	OPTIKA PRO VIEW	OPTIKA LITE VIEW
Simultaneous management of seve	Simultaneous management of several cameras		Х
IMAGE acquisition		X	X
	tiff	х	Х
	jpg	х	Х
	bmp	X	X
Image formats	png	X	X
	рсх	х	X
	jp2	x	X
	dcm	x	X
IMAGE acquisition		x	X
	avi	x	X
	wmv	X	X
	mp4	X	X
VIDEO formats	asf	x	X
	3gp	X	X
	mov	X	X
	h264	X	X
Continuous automatic exposure Manual Exposure	h265	X	X
Continuous automatic exposure	Continuous automatic exposure		Х
	Manual Exposure		Х
Mobile spot for exposure Resizable spot for exposure Colour acquisition	Mobile spot for exposure		X
Resizable spot for exposure	Resizable spot for exposure		X
Colour acquisition	Colour acquisition		X
Grey-scale acquisition	Grey-scale acquisition		X
Manual Time-Lapse	Manual Time-Lapse		
Automatic Time-Lapse	Automatic Time-Lapse		
Fast Image Acquisition	Fast Image Acquisition		X
Focus Indicator	Focus Indicator		
White Balance	White Balance		X
Black balance	Black balance		
Background correction	Background correction		
Dark Field Correction	Dark Field Correction		X
Image Enhancement	Image Enhancement		X
Live Histogram		Х	Х
Elio	Horizontal	Х	X
Flip	Vertical	Х	Х
Rotate		х	

	FUNCTION	OPTIKA PRO VIEW	OPTIKA LITE VIEW
PROCESSING	Several function of image processing (filters)	x	
	Multiple image combining	x	
	EDF (Extended Depth of Focus)	×	
	Colour Combine (Multi-Fluorescence Imaging)	×	
	Shift Correction	×	
	HDR (High Dynamic Range)	×	
	Layer Management	x	
	Text Overlay	x	
	Ruler Overlay	×	
	Measurement Overlay	×	
	Grids	×	

4083.Wifi, **4083.4** and **4083.EC2** work with Vision Lite only.

C-E2 work with PROVIEW and LITEVIEW only.

Cameras with HDMI connection only, do not require any software; C-HP works with PROVIEW only, isn't compatible with Linux and MAC OS.